

# Application Sheet

## RUBELLA IgG

**REF 51208**

## Instrument: Elisys Uno

**REF 17350**

### Setting:

**RUBG\_51208\_R6**

For all essential product information, please refer to the User Manual of the analyzer and the instructions for use for reagents, controls, calibrators and the Certificate of Analysis. Please refer to the instruction for use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyzer.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

### Material Required

Material	REF	Comments
Rubella IgG	51208	
Reagent Rack, Infect. Diseases & Steroids & Autoimmune Elisys Uno	17359/2	all reagents
Sample Rack (12 mm) Elisys Uno or Sample Rack (13 mm) Elisys Uno	17359/7 or 17359/8	Sample and dilution; according to sample collection tube
Sample vials 2 ml	17350/16	Sample dilution

### Additional Notes

Although the results should be interpreted with regard to a cut-off value, the assay is programmed as a quantitative assay to get concentrations for positive samples in the report. For this reason all results are shown with concentrations. Concentrations below 10 U/ml and above 65 U/ml are calculated by extrapolation and should be considered as "out of range".

Due to the usable space on the racks the number of samples that can be processed is limited to 48.