

# Application Sheet

## HSV 2 IgG

REF 51226

## Instrument: Elisys Uno

REF 17350

### Setting:

HSV2G\_51226\_R4

For all essential product information, please refer to the User Manual of the analyser and the Instructions for Use for reagents, controls and calibrators. Please refer to the Instructions for Use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyser.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

### Material Required

| Material   | REF                | Comments   |
|--|--------------------|--|
| HSV 2 IgG  | 51226              |  |
| Reagent Rack, Infect. Diseases & Steroids                              | 17359/2            | all reagents   |
| Sample Rack (13 mm) Elisys Uno<br>or<br>Sample Rack (12 mm) Elisys Uno | 17359/8<br>17359/7 | Sample and dilution, according to sample collection tube |
| Sample vials 2 ml  | 17350/16           | Sample dilution  |

### Additional Notes

For technical reasons the assay step „start incubation time” requires to define a corresponding reagent. This reagent only exists in the software, there is no reagent bottle for it. However, this virtual reagent needs to be allocated to a position on the reagent rack, any position can be chosen. Do not place any other reagent bottle on this position while the virtual bottle is allocated to it.

Due to the usable space on the racks the number of samples that can be processed is limited to 48.