

# Application Sheet

## HBsAg ultra sens (96T) Human

REF 51050

## Instrument: HumaReader Single plus

REF 18000

### Setting:

HBsAg\_51050\_R2

For all essential product information, please refer to the User Manual of the analyzer and the instructions for use for reagents, controls, calibrators and the Certificate of Analysis. Please refer to the instruction for use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyzer.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

### Additional Notes

The reading may be cancelled automatically because of an invalid NC-value. In this case read the absorbance values only (without automated interpretation). If two NC-values meet the validity criterion use these values to calculate the cut-off and interpret the test manually. The criteria for LPC must always be checked by the user. The criterion for PC is checked automatically.

HumaReader Single plus requests a sample order that differs from the IFU. NC in wells A1-C1, PC in well D1.