

Application Sheet

Toxo IgM μ -Capture

REF 51119

Instrument: HumaReader Single plus

REF 18000

Setting:

TOXOMu 51119 R4

For all essential product information, please refer to the User Manual of the analyser and the Instructions for Use for reagents, controls and calibrators. Please refer to the Instruction for Use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyser.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

Additional Notes

Please note that the reader requires strict positioning of controls. PCH cannot be recognized by the reader as an additional positive control. For this reason, PCH has to be loaded as a sample.

Well Position	Layout of test
A1/B1	Blank
C1/D1	NC
E1/F1	PCL as 'PC'
G1/H1	PCH as 'sample 1'
A2...	samples

The validity criteria "MPCH > MPCL" and "MPCL/COV \geq 3.0" must be verified manually.