

Application Sheet

fPSA

REF 52035

Instrument: HumaReader HS

REF 16670

Setting:

fPSA_52035_R2

For all essential product information, please refer to the User Manual of the analyzer and the instructions for use for reagents, controls, calibrators and the Certificate of Analysis. Please refer to the instructions for use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyzer.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

Additional Notes

The validation criterion $CAL F \geq 1.3$ is checked automatically. If the criterion fails, no result will be displayed.

The criterion "The difference between the duplicates of CAL F does not exceed 10%" is not checked by the reader.