

Application Sheet

Toxo IgM μ -Capture

REF 51119

Instrument: HumaReader HS

REF 16670

Setting:

TOXO_M_ μ C_51119_R4

For all essential product information, please refer to the User Manual of the analyser and the Instructions for Use for reagents, controls and calibrators. Please refer to the Instruction for Use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyser.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

Additional Notes

PCL has to be defined as 'PC1' and PCH as 'PC2'.

The validation criteria are checked automatically. If the criteria fail, no result will be displayed.