

# HumaTurb Solution

## Molecular Screening of *MTBC* and other Pathogens in Regional and Reference Laboratories

- Adaptable set up for different throughputs
- Multiple parameters can be tested simultaneously
- Flexible data transfer via USB and LIS connectivity




# HumaTurb Solution

## Incubator for Amplification & Measurement of Turbidity in Real-time

### Features and technical specifications

- HumaTurb C – Control unit with operation panel and built-in printer
- HumaTurb A – Amplification unit with 2 independent operation reaction blocks
- HumaHeat – Incubator for the sample lysis of the Loopamp™ PURE heating tubes
- HumaTurb C is connectable with up to 6 HumaTurb A units for increased throughput
- Flexible solution for up to 96 tests/run
- Automated result reading



981250/2019-03 © 2019 HUMAN

**Human**

Diagnostics Worldwide