

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml

Version 02/ 09-2024

371 Level 1 12x5ml

Level 1 40441

REF

372 Level 2 12x5ml

LOT

40440

Level 2 40442

IVD



2026-02-28

373 Level 3 12x5ml

Level 3 40443

370X MiniPak 3x5ml



If a reagent kit with a new LOT number is used, please check the link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>

REF	Assay Parameter Ensayo Essai	Valid from Gültig ab Válido a partir de Valable à partir de RGT Lot	Units Einheiten Unidades Unité	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle
84600	TSH	23001	µIU/ml	0.530	0.345 - 0.716	5.90	3.84 - 7.97	39.6	25.8 - 53.5
84620	fT3	23002	pg/ml	2.93	1.90 - 3.96	7.25	4.71 - 9.79	10.5	6.85 - 14.2
84615	T3	23001	ng/ml	0.930	0.605 - 1.26	2.42	1.57 - 3.27	3.66	2.38 - 4.94
84610	fT4	23001	pg/ml	8.66	5.63 - 11.70	20.7	13.50 - 28.0	33.6	21.9 - 45.4
84605	T4	23001	ng/ml	57.2	37.2 - 77.3	117	76.1 - 158	157	102 - 212
82120	FSH	23001	mIU/ml	7.27	4.73 - 9.81	23.4	15.2 - 31.5	51.3	33.4 - 69.3
82115	LH	23001	mIU/ml	2.02	1.31 - 2.73	26.2	17.1 - 35.4	89.5	58.2 - 121
83130	PSA	23001	ng/ml	0.790	0.514 - 1.07	3.69	2.40 - 4.98	25.0	16.2 - 33.7
83125	fPSA	23001	ng/ml	0.520	0.338 - 0.702	1.62	1.05 - 2.19	9.46	6.15 - 12.8
83100	AFP	23001	ng/ml	19.8	11.9 - 27.7	105	63.0 - 147	303	182 - 425
82100	HCG	23001	mIU/ml	4.46	2.90 - 6.02	26.7	17.4 - 36.0	193	125 - 261
82105	Progesterone	23001	ng/ml	0.150	<0,100 - 0.203	12.2	7.94 - 16.5	29.2	19.0 - 39.5
82110	Testosterone	23001	ng/ml	0.720	0.432 - 1.01	3.4	2.04 - 4.76	6.80	4.08 - 9.52
85000	25-OH Vitamin D	23001	ng/ml	20.2	13.1 - 27.3	26.3	17.1 - 35.5	21.0	13.6 - 28.3
83105	CEA	23001	ng/ml	2.15	1.40 - 2.90	16.6	10.8 - 22.4	32.9	21.4 - 44.4
82130	Estradiol	23001	pg/ml	49.5	32.2 - 66.8	208	135 - 281	694	451 - 936
85049	Ferritin*	23001	ng/ml	76.5	49.7 - 103	188	122 - 254	425	277 - 574
82125	Prolactin	23001	ng/ml	9.60	6.24 - 13.0	27.3	17.8 - 36.9	57.9	37.6 - 78.1
85055	Vitamin B12	23001	pg/ml	441	287 - 596	673	438 - 909	944	614 - 1275

* No CE - only for use outside Europe