

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

REF 370, **LOT** 40470 Level 1 - 3  2028-01-31

Revision 1

REF		Units Einheiten Unitades Unitées	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Rango Marge
HUMAN Hormones & Tumor Markers								
			Level 1 - 40471		Level 2 - 40472		Level 3 - 40473	
52010	AFP	ng/ml	18.7	12.5 - 24.9	93.8	62.9 - 125	241	162 - 322
52020	CEA	ng/ml	2.04	1.36 - 2.71	12.0	8.1 - 16.0	20.0	13.4 - 26.7
55050	Cortisol	ng/ml	37.3	23.9 - 50.7	208	133 - 283	365	233 - 496
55060	DHEA-S	µg/ml	0.61	0.33 - 0.88	1.47	0.81 - 2.13	4.32	2.38 - 6.27
55030	Estradiol	pg/ml	167	93.5 - 240	454	254 - 654	874	4.89 - 1258
55040	Estriol	ng/ml	2.70	1.49 - 3.91	7.70	4.23 - 11.2	14.4	7.91 - 20.8
52035	fPSA	ng/ml	0.57	0.38 - 0.76	2.92	1.95 - 3.88	> 10	10.0
53020	FSH	IU/l	6.78	4.54 - 9.01	22.7	15.2 - 30.1	46.3	31.0 - 61.5
54015	ft3	pg/ml	2.65	1.78 - 3.52	5.36	3.59 - 7.13	7.90	5.29 - 10.5
54025	ft4	ng/dl	0.75	0.50 - 1.00	2.01	1.35 - 2.67	3.27	2.19 - 4.35
53040	hCG	IU/l	5.94	3.98 - 7.90	19.9	13.4 - 26.5	223	150 - > 250
53010	LH	IU/l	1.99	1.33 - 2.64	23.4	15.7 - 31.1	79.9	53.5 - 106
53030	Progesterone	ng/ml	0.94	0.52 - 1.36	10.4	5.74 - 15.1	31.8	17.5 - > 40
55020	Prolactin	ng/ml	8.72	5.85 - 11.60	23.5	15.8 - 31.3	57.4	38.5 - 76.4
52030	PSA	ng/ml	2.35	1.29 - 3.40	8.00	5.20 - 10.8	26.2	17.1 - 35.4
54010	T3	ng/ml	0.98	0.66 - 1.30	2.40	1.61 - 3.20	4.06	2.72 - 5.40
54020	T4	µg/dl	5.09	3.41 - 6.76	9.13	6.12 - 12.1	16.0	10.8 - 21.3
55010	Testosterone	ng/ml	1.38	0.8 - 1.97	3.62	2.10 - 5.14	7.38	4.28 - 10.5
54030	TSH	mIU/l	0.182	0.10 - 0.269	4.86	3.16 - 6.56	30.9	20.1 - 41.7

* The marked level is not applicable for the indicated parameter as the target value is out of the measuring range of the assay.



Diagnostics Worldwide

Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany

Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de