

HUMAN

Target Values

BioRad Liquichek Cardiac Markers Plus Control

180 Trilevel 6×3ml

Version 01/ 03 - 2022

181 Level 1 6×3ml

Level 1 87831

REF

182 Level 2 6×3ml

LOT

87830

Level 2 87832

IVD



2024-03-31

183 Level 3 6×3ml

Level 3 87833

180X MiniPak 3×3ml

Link for latest version of target value sheet:

<https://www.human.de/cli/target-value-sheets/>



REF	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de	Units Einheiten Unidades	Target Sollwert Valor asignado Cible	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	Range Bereich Intervalo
	Essai	Valable à partir de RGT Lot	Unité	Intervalle	Intervalle	Intervalle	Intervalle	Intervalle	Intervalle
				Level 1 - 87831		Level 2 - 87832		Level 3 - 87833	
85030	Troponin I	21001	ng/ml	0.55	0.41 - 0.69	4.30	3.23 - 5.38	16.34	12.26 - 20.43
85040	CK-MB	21001	ng/ml	3.61	2.71 - 4.51	15.93	11.95 - 19.91	79.15	59.36 - 98.94
85035	Myoglobin	21001	ng/ml	62.31	46.73 - 77.89	193.88	145.41 - 242.35	464.30	348.23 - 580.38

#87830
Rev. 01 | valid of 24.03.2022



Human Gesellschaft für Biochemia und Diagnostica mbH
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de