

HUMAN

Target Values

BioRad Liquichek Cardiac Markers Plus Control

180 Trilevel 6×3ml

Version 01/ 07 - 2024

181 Level 1 6×3ml

Level 1 87871

REF

182 Level 2 6×3ml

LOT

87870

Level 2 87872

IVD

🕒 2026-04-30

183 Level 3 6×3ml

Level 3 87873

180X MiniPak 3×3ml

If a reagent kit with a new LOT number is used, please check the link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>



REF	Assay Parameter Ensayo Essai	Valid from Gültig ab Válido a partir de Valable à partir de [RGT] Lot	Units Einheiten Unidades Unité	Target		Range Bereich Intervalo Intervalle	Target		Range Bereich Intervalo Intervalle	Target	
				Sollwert	Valor asignado Cible		Sollwert	Valor asignado Cible		Sollwert	Valor asignado Cible
				Level 1 - 87871		Level 2 - 87872		Level 3 - 87873			
85030	Troponin I	24001	ng/ml	0.16	0.12 - 0.20	1.04	0.78 - 1.30	3.61	2.71 - 4.51		
85040	CK-MB	23002	ng/ml	3.22	2.42 - 4.03	13.41	10.06 - 16.76	66.39	49.79 - 82.99		
85035	Myoglobin	23002	ng/ml	59.61	44.71 - 74.51	202.91	152.18 - 253.64	490.72	368.04 - 613.40		

