


HbA1c liquidirect Controls

LOT 25001

Version 1 /05 - 2024

REF 10775

 2027-02-28

Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

Device Gerät Instrumento Equipement		Unit Einheit Unidad Unité	Control Blood normal CBN				Control Blood abnormal CBA			
			Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs			Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs		
HumaStar 300SR	mmol/mol Hb		35.2	28.7	—	41.7	109	98.0	—	121
	%		5.4	4.8	—	6.0	12.1	11.1	—	13.2
HumaStar 600	mmol/mol Hb		43.4	36.6	—	50.3	116	105	—	127
	%		6.1	5.5	—	6.8	12.8	11.8	—	13.8
other	mmol/mol Hb		35.2	28.7	—	41.7	109	98.0	—	121
	%		5.4	4.8	—	6.0	12.1	11.1	—	13.2

CS-HbA1c

Rev. 007 | valid of 15.11.2022

INF 1077502 - Version 1 /05 - 2024



Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 · 65205 Wiesbaden · Germany
Telefon +49 6122-9988-0 · Telefax +49 6122-9988-100 · e-Mail human@human.de