

TURBIDOS

REF 13010

LOT

24001

Version 1 /08 - 2024



2026-04-30

Temporary target values of RF valid until

2024-12-31

Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

	Method Methode Método Méthode	SI Unit SI Einheit Unidad SI Unité SI	Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs	Unit Einheit Unidad Unité	Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs
HUMAN reagent kits							
Level 1							
RHEUMATOID FACTORS from LOT 23501							
HumaLyzer photometer HumaStar analyzer	Immunoturbidimetry	kIU/l	43.9	31.4 - 57.0	IU/ml	43.9	31.4 - 57.0
Level 2							
RHEUMATOID FACTORS from LOT 23501							
HumaLyzer photometer HumaStar analyzer	Immunoturbidimetry	kIU/l	98.7	76.2 - 122	IU/ml	98.7	76.2 - 122