

# HUMAN

Target Values

## BioRad Lyphocheck Immunoassay Plus Control

**REF** 370, **LOT** 40440 Level 1 - 3 2026-02-28

Revision 1

REF		Units Einheiten Unitades Unitées	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Rango Marge
<b>HUMAN Hormones &amp; Tumor Markers</b>								
			Level 1 - 40441		Level 2 - 40442		Level 3 - 40443	
52010	AFP	ng/ml	13.7	9.20 - 18.3	87.2	58.4 - 116	282	189 - 375
52020	CEA	ng/ml	1.07	0.71 - 1.42	11.4	7.6 - 15.1	18.2	12.2 - 24.2
55050	Cortisol	ng/ml	37.6	24.1 - 51.2	213	136 - 289	341	218 - 464
55060	DHEA-S	µg/ml	0.60	0.33 - 0.87	1.78	0.98 - 2.58	4.56	2.51 - 6.62
55030	Estradiol	pg/ml	196	110 - 282	530	297 - 764	1206	676 - 1737
55040	Estriol	ng/ml	2.01	1.11 - 2.92	6.27	3.45 - 9.09	15.1	8.30 - 21.9
52035	fPSA	ng/ml	< 2,00	-	4.33	2.9 - 5.76	26.8	17.9 - 35.6
53020	FSH	mIU/ml	6.79	4.55 - 9.03	22.1	14.8 - 29.3	50.5	33.8 - 67.1
54015	fT3	pg/ml	2.77	1.86 - 3.69	5.38	3.61 - 7.16	7.54	5.05 - 10.0
54025	fT4	ng/dl	0.92	0.62 - 1.23	2.07	1.39 - 2.76	2.72	1.82 - 3.61
53040	hCG	mIU/ml	6.29	4.21 - 8.36	21.8	14.6 - 29.0	173	116 - 230
53010	LH	mIU/ml	1.79	1.20 - 2.38	22.8	15.3 - 30.4	65.4	43.8 - 86.9
53030	Progesterone	ng/ml	0.91	0.50 - 1.32	12.2	6.72 - 17.7	39.3	21.6 - > 40
55020	Prolactin	ng/ml	6.90	4.62 - 9.18	18.3	12.2 - 24.3	47.1	31.5 - 62.6
52030	PSA	ng/ml	1.84	1.01 - 2.67	5.70	3.70 - 7.69	25.1	16.3 - 33.9
54010	T3	ng/ml	1.00	0.67 - 1.33	2.56	1.71 - 3.40	3.81	2.55 - 5.06
54020	T4	µg/dl	5.25	3.52 - 6.99	8.82	5.91 - 11.7	12.4	8.32 - 16.5
55010	Testosterone	ng/ml	1.44	0.83 - 2.04	4.70	2.72 - 6.67	9.60	5.57 - 13.6
54030	TSH	mIU/l	0.21	0.12 - 0.31	3.70	2.40 - 4.99	25.6	16.6 - 34.6

\* The marked level is not applicable for the indicated parameter as the target value is out of the measuring range of the assay.

