

HUMAN

Target Values

BioRad Lyphocheck Tumor Marker Plus Control

REF 368X, **LOT** 74620 Level 1 - 3 ⌚ 2025-05-31

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten	Sollwert	Bereich	Sollwert	Bereich	Sollwert	Bereich
		Unités	Valor meta	Rango	Valor meta	Rango	Valor meta	Rango
			Valeur souhaitée	Marge	Valeur souhaitée	Marge	Valeur souhaitée	Marge
HUMAN Tumor Markers								
			Level 1 - 74621		Level 2 - 74622		Level 3 - 74623	
52010	AFP	ng/ml	7.53	5.04 - 10.0	77.2	51.7 - 103	249	154 - 345
52020	CEA	ng/ml	3.35	2.24 - 4.46	28.1	18.8 - 37.3	100.0	67.3 - 134
52050	CA125	U/ml	19.7	13.2 - 26.2	59.1	39.6 - 78.6	205	137 - 273
52080	CA15-3	U/ml	33.6	20.8 - 46.3	78.7	52.7 - 105	157	105 - 209
52060	CA19-9	U/ml	20.5	13.70 - 27.2	41.4	27.7 - 55.1	103.0	68.9 - 137
52035	fPSA	ng/ml	< 0,3		5.38	3.60 - 7.16	> 10	
53040	PSA	ng/ml	0.22	0.12 - 0.32	10.4	6.78 - 14.1	34.0	22.1 - 45.9

