

HUMAN

Target Values

BioRad Lyphocheck Tumor Marker Plus Control

REF 368X, **LOT** 74650 Level 1 - 3 ⌚ 2026-10-31

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten	Sollwert	Bereich	Sollwert	Bereich	Sollwert	Bereich
		Unités	Valor meta	Rango	Valor meta	Rango	Valor meta	Rango
			Valeur souhaitée	Marge	Valeur souhaitée	Marge	Valeur souhaitée	Marge
HUMAN Tumor Markers								
			Level 1 - 74641		Level 2 - 74642		Level 3 - 74643	
52010	AFP	ng/ml	6.44	4.32 - 8.57	78.0	52.2 - 104	238	160 - 317
52020	CEA	ng/ml	2.89	1.94 - 3.85	24.2	16.2 - 32.2	81.3	54.4 - 108
52050	CA125	U/ml	19.8	13.2 - 26.3	59.3	39.7 - 78.9	151	101 - 201
52080	CA15-3	U/ml	23.6	15.8 - 31.4	50.5	33.9 - 67.2	110	73.5 - 146
52060	CA19-9	U/ml	12.4	8.28 - 16.4	27.5	18.4 - 36.5	72.6	48.6 - 96.5
52035	fPSA	ng/ml	< 0,3		2.89	1.94 - 3.85	> 10	
53040	PSA	ng/ml	0.260	0.145 - 0.382	10.93	7.11 - 14.8	31.3	20.4 - 42.3

