

**HUMAN**  
Target Values  
**BioRad Lyphocheck Tumor Marker Plus Control**  
**REF 368X, LOT 74660 Level 1 - 3 2027-07-31**

REF		Units Einheiten Unitades Unitées	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
<b>HUMAN Tumor Markers</b>								
			Level 1 - 74661		Level 2 - 74662		Level 3 - 74663	
52010	AFP	ng/ml	8.00	5.36 - 10.6	69.8	46.8 - 92.9	221	148 - 294
52020	CEA	ng/ml	2.69	1.80 - 3.58	21.9	14.7 - 29.1	81.8	54.8 - 109
52050	CA125	U/ml	18.9	12.7 - 25.1	56.5	37.8 - 75.1	164	110 - 218
52080	CA15-3	U/ml	30.7	20.5 - 40.8	56.5	37.9 - 75.2	127	85.3 - 169
52060	CA19-9	U/ml	15.2	10.20 - 20.2	34.0	22.8 - 45.2	90.4	60.6 - 120
52035	fPSA	ng/ml	< 0,3		2.94	1.97 - 3.92	> 10,0	
53040	PSA	ng/ml	0.20	0.11 - 0.29	12.3	7.98 - 16.6	34.5	22.4 - 46.6

