

Standard RHEUMATOID FACTORS

Version 1/02 - 2024

REF

11361

LOT

24002



2025 - 07 - 31

Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

Calibration Values / Kalibrationswerte / Valores de calibración / Valeurs de calibration

Device Gerät Instrumento Equipement	Reference material/method Referenzmaterial/-methode Material/método de referencia Matériel/méthode de référence	Unit Einheit Unidad Unité	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée
HUMAN RHEUMATOID FACTORS - Immunturbidimetry Latex enhanced			
Multipurpose reagent on HumaStar 100 & 200 LOT ≥ 23501	In house Master calibrator traceable to 1 st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
Multipurpose reagent on HumaLyzer 4000 LOT ≥ 23501	In house Master calibrator traceable to 1 st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
Multipurpose reagent for other LOT ≥ 23501	In house Master calibrator traceable to 1 st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
HumaStar System reagent on HumaStar 300SR LOT ≥ 23501	In house Master calibrator traceable to 1 st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
HumaStar System reagent on HumaStar 600 LOT ≥ 23501	In house Master calibrator traceable to 1 st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180

CA_RF_11361

Rev. 002 | valid of 15.11.2022

INF 1136102 - Version 1/02 - 2024



Human

Human Gesellschaft für Biochemica und Diagnostica mbH

Max-Planck-Ring 21 · 65205 Wiesbaden · Germany

Telefon +49 6122-9988-0 · Telefax +49 612-9988-100 · e-Mail human@human.de